# INSIZE

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



### 4. KEY DIAGRAM


## 5. IMPORTANT GUIDELINES

- » Place the moisture analyzer on a stable, horizontal surface.
- » Avoid placing near radiators or direct sunlight to prevent heat interference.
- » Protect from direct airflow (open windows, doors).
- » Avoid vibrations during weighing.
- » Keep the instrument away from high humidity, steam, and dust.
- » If moving the device from cold to warm environments, allow it to adjust for 2 hours before use to prevent condensation.
- » Static electricity and electromagnetic fields can affect accuracy. Ensure stable power and adjust location if needed.

## 6. OPERATION INSTRUCTIONS

### Adjust Level:

Align the instrument using the leveling bubble until it is within the designated circle.

Note: Regularly check whether the instrument is in a level.

Power On: Turn on the device and allow 20-30 minutes for preheating.

### Calibration:

#### 8701-RM20

Hold [Tare] button until "LOAD 100" appears on the display, place a 100 g weight on the scale, and hold [Tare] button again for about 3-5 seconds. Calibration completes when "100.000 g" appears on the display. After successful calibration, the instrument automatically returns to the weighing mode.

#### 8701-RM30 / 8701-RM40

Press [Tare] button to clear, then press [Cal] button, "LOAD 100" appears on the display. Take off the scale and place a 100 g weight on the triangle support frame (remove sample pan), and hold [Cal] button for about 3-5 seconds. Calibration completes when "100.000 g" appears on the display. After successful calibration, the instrument automatically returns to the weighing mode. Place back the scale.

Zeroing: Press [Tare] button and wait until the display shows zero.

Sample Weighing: Press [Tare] button to zero the scale, then place the sample, and wait for the "OK" signal on the lower left side of the display to ensure stability. After that the displayed value on the display is the sample weight value.

Overload Warning: Exceeding the maximum load will display an error ("E" or "H"). Reduce the load to prevent damage.

Settings: Press [Menu] button to enter settings. Navigate with [Tare] button, confirm with [Start] button, and exit by pressing [Tare] button until "ESC" appears, then confirm with [Start] button.

## Menu

- Factory setting\*      - User setting
  
- STOP      - CLOC (Time shutdown)\*
  - AUTO (Automatic weightless shutdown)
  - ESC
  
- LOSE Weight-loss rate setting      - A1
  - A2\*
  - A3
  - A4
  - ESC

Note: Those marked with \* are factory defaults

## 7. FIRST TEST

In the first measurement, the equipped sample (glass fiber filter paper with water) is used to measure the water content. During the first measurement, the operation is according to factory settings.

- » Connect the moisture analyzer to the power supply and wait until "0.000 g" displays.
- » Press [Tare] button, open the heating chamber, place the sample, close the heating chamber, and press [Start] button to do a drying test according to the default settings. You can monitor the process on the display (temperature, drying time, weight).
- » Display the current temperature of the heating unit, the past drying time, and the current mass or moisture content after drying (press the unit conversion button [Unit] to freely switch between mass and % content).
- » Press [Stop] button to end the test at any time. The final moisture content will be displayed.
- » Press the [Tare] button to return to the normal weighing state and remove the sample.

! Burn hazard: Samples, sample trays, and brackets are now very hot! After the drying process is completed, open the heating chamber and gently remove the sample tray using the sample tray bracket.

## 8. TESTING PROCESS

- » Set heating method, time, temperature, and shutdown mode.
- » Prepare the sample:
  1. Place the empty sample tray on the sample tray bracket
  2. Press the [Tare] button to remove the weight of the sample tray
  3. Take out the sample tray and place the sample, preferably greater than 0.5 g
  4. Distribute the sample evenly on the sample tray
- » Place the sample tray containing the sample back on the bracket, and the display screen displays the weight of the sample. Close the heating chamber.
- » Press the [Start] button to start the measurement.
- » The display screen displays the current weight value, and press the [Unit] button to convert the current measurement result.
- » After the measurement is completed, the display screen displays the final sample moisture content, press the [Tare] button to return to the normal weighing mode.
- » Press the [Print] key to send the current measurement results to a printer or other peripheral device.

Note: You can press the [Stop] button at any time during the test to terminate this measurement.

## Shutdown Mode Settings

The shutdown mode is used to determine the automatic test end point of the moisture analyzer and display the results. The moisture analyzer has two shutdown modes to set: timed shutdown and automatic shutdown.

» Timed shutdown: manually set the test time.

Press the time button and the time symbol flashes. The hydrometer displays the default time of the system or the time set last time (minutes).

Reduce the temperature setting for 1 minute every time you press Tare, increase the temperature setting for 1 minute every time you press Start. Press [Time] button to confirm the save.

» Automatic shutdown: automatically set the test time.

Press [Menu] to display AUTO.

Press the "Start" button to confirm that the setting is successful when the letter "OK" appears in the lower left corner of the display, and the instrument will return to the weighing state.

(If you want to change the automatic mode to manual mode, press [Menu] button the instrument displays "AUTO", then press the [Tare] button to find "CLOCK" and press the [Start] button to confirm, the setting is successful when "OK" appears in the lower left corner of the display, and the instrument will return to weighing mode.)

## Loss Rate Settings

» A1 (5mg/20s): Rapid trend detection.

» A2 (5mg/50s): Quick drying.

» A3 (5mg/60s): Standard setting for most samples.

» A4 (5mg/90s): For medium-dry samples.

» Press [Menu] button to enter the option and press the [Tare] button to select.

» Press the [Start] button to confirm, and the instrument will display the graphic symbol of automatic shutdown and return to the normal weighing state.

## 9. WARNING

The following substances may pose a risk of fire, explosion, damage, or injury. For substances with potential safety hazards, please carefully analyze the possible hazardous consequences. In this case, the instrument should be watched by a dedicated person and wear corresponding protective devices.

### 1. Volatile substances

For volatile substances, it is recommended to use a rapid heating method to limit the loss of water before the sample evaporates. The test should be conducted in a safe and dry environment with a sufficiently low temperature to prevent explosion or fire. Please use 1 g or less of sample for testing.

### 2. Toxic substances

Contains toxic or corrosive substances and should be dried in ventilated shop windows.

### 3. Corrosive substances

Heating substances that produce corrosive gases, such as acidic substances, should be tested with as few samples as possible. Because of its production.

## 10. MAINTENANCE


1. Disconnect the power supply when cleaning the moisture analyzer.
2. Ensure no liquid enters the moisture analyzer.
3. Please ensure that the moisture analyzer is in the cooling state before cleaning.
  - » Clean the instrument regularly.
  - » The cover surface and temperature sensor can be slightly dipped in a lint-free cloth or neutral cleaning agent scrub.
  - » Glass surface can be cleaned using glass cleaning agent.
  - » Do not use solvents, irritating chemicals such as ammonia and abrasive solvents.

## 11. RS232C

This instrument is equipped with RS232C interface for connecting peripherals (printers, computers).

In order to successfully connect the instrument to peripheral equipment, the following conditions must be met:

- » Connect instruments using cables that match peripheral interfaces.
- » Communication parameters (baud rate, data bits and parity bits) of the instrument and peripheral devices must match. This data exchange is asynchronous using ASCII codes.

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**INSIZE is a global manufacturer of measuring technology** bringing technological innovations with representation in 75 countries around the world. Measuring instruments of the INSIZE brand represent the optimal solution without compromises to meet even the most demanding measuring needs.

**INSIZE measuring instruments will pleasantly surprise you with:**  
| high-quality craftsmanship | reliable performance | affordable prices

INSIZE provides a complete portfolio of measuring instruments » calipers, height gauges, dial indicators, micrometers, roughness testers, hardness testers, measuring microscopes, optical measuring devices, profile projectors, testing machines, contour gauges, roundness measuring machines, thickness gauges, gages, protractors, force gauges, meters, scales, videoscopes, torque wrenches and accessories including stands for measuring instruments.

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