INSIZE

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1. FEATURES

- » High accuracy and sensitivity.
- » Automatic correction function.
- » Unit conversion function allows switching between “g” (grams), “ct” (carats), and “ozt” (ounces).
- » The scale is equipped with a data transfer interface, which can be connected to a printer for data printing or to a computer for data collection and analysis.
- » The counting function helps users easily count large quantities of objects.


2. TECHNICAL SPECIFICATION

Code	8301-300	8301-600	8301-1200	8301-3000	8301-6000
Maximum weighing capacity	300 g	600 g	1200 g	3000 g	6000 g
Minimum weighing capacity	0,2 g	0,2 g	0,2 g	2 g	2 g
Resolution (d)	0,01 g	0,01 g	0,01 g	0,1 g	0,1 g
Verification interval (e)	0,1 g	0,1 g	0,1 g	1 g	1 g
Accuracy	m ≤ 50 g: ±0,05 g 50 g < m ≤ 200 g: ±0,1 g m > 200 g: ±0,15 g			m ≤ 500 g: ±0,5 g 500 g < m ≤ 2000 g: ±1 g m > 2000 g: ±1,5 g	
Weighing platform size	∅135 mm			190×160 mm	
Operating temperature	0-40 °C				
Operating humidity	≤ 80 % RH				
Power supply	220 V, 50/60 Hz				
Dimensions (L×W×H)	215×196×62 mm				

Note: Models 8301-300, 8301-600, and 8301-1200 are supplied with a calibration weight and protective glass.

3. STRUCTURE

1. Display
2. Combined button
3. Calibration button
4. TARE button
5. Counting button
6. Power switch
7. Weighing platform
8. Power socket
9. Fuse holder
10. Serial port RS232
11. Leveling bubble

Accessories:

Protective glass and calibration weight for models 8301-300 to 8301-1200.

4. OPERATION

1. Turn on the power, set the plane, and switch on the scale. The display will show "F----1" to "F----9," then "0." Ideally, the scale should be switched on for at least 15 minutes before use to allow it to warm up. It is normal for the scale to show slight deviations when first turned on, which will stabilize after a while.

2. TARE button: If the display shows a deviation from zero without any load, press the "TARE" button to reset the scale to zero. To remove the weight of a container, place the container on the weighing platform, wait for the value to stabilize, and then press the "TARE" button. The scale will display zero weight. You can then place the object to be weighed into the container. The weight displayed on the screen at this time is the net weight of the object in the container.

3. "CAL" Calibration Button: If the scale is newly purchased or has not been used for a long time, it must be calibrated. After warming up the scale (minimum 15 minutes), press the "CAL" button. The display will show "C XXX" (where XXX is the weight of the calibration mass, e.g., "C 200"). Place the calibration weight, and after stabilization, the display will show the weight and the stabilization symbol "g." The calibration is complete at this point. If "C----F" appears, it means the zero point is unstable—press the "TARE" button and repeat the calibration. If "F----H" appears, it means the calibration weight exceeds the scale's maximum capacity. Repeat the calibration with a weight within the scale's range.

4. Using the "COUNT" button:

1) Selecting object quantity: For accurate counting, first select the number of objects (1-10-20-50-100). For very light objects and those with very similar weights, select the largest number of objects.

2) If the scale is empty, place one or more objects on the scale, and it will display the objects' weight. Press the "COUNT" button, and the scale will display "1," with the "pcs" indicator lighting up on the right side of the display, indicating that the scale is in counting mode. The inserted objects are counted as 1 unit. If the inserted objects are not 1, press the "COUNT" button again to select the number corresponding to the inserted objects. Then place the same objects, and the displayed value is the total number of objects. Pressing the "COUNT" button again returns the scale to weighing mode.

5. Combined button (unit conversion and printing):

For detailed settings, refer to the section "Setting method and meaning of the fifth point of the electronic scale." The default factory setting is the unit conversion button.

1) In weighing mode, press this button to switch between three weight units: "g" (gram), "ct" (carat), and "ozt" (ounce). The selected unit is displayed in the lower right corner of the screen.

2) You can also press this button with an object already on the platform, and the units will convert.

6. Data output: The scale is equipped with an RS232C interface that can connect to a printer or computer. For computer communication, the software must have the following settings:

1) Communication protocol: Baud rate: adjustable between 600 and 9600, factory default is 600; data bits: 8 bits; stop bit: 1 bit; no parity bit.

2) Output data: Press the print button once (see Appendix C5 for print button settings), or the serial port will receive a print command for data output once. Data is transmitted in a 14-bit ASCII format.

symbol	data	unit	enter	wrap
1 bit	8 bit	3 bit	1 bit	1 bit

If the data is not negative, the first position is a space, and the unshown data is represented as a space.
Example: If 123.45 g is displayed, the output data for this measurement will be `□□□123.45□g□↓←`. If -1234.5 g is displayed, the output data for this measurement will be `-□□1234.5□g□↓←`. Total 14-bit data.

3) Receive command: First, the corresponding signal must be received, and then the command is accepted. When signal 27 (factory state), which is 1B in hexadecimal, is passed, the following command is sent (data is in hexadecimal):

- » 1B 70: Print (requires the scale to send data once)
- » 1B 71: Calibration
- » 1B 72: Counting
- » 1B 73: Unit conversion
- » 1B 74: Tare

4) Pin definitions of RS232C output interface (DB9 socket (female type)): Pin 2: RXD Pin 3: TXD Pin 5: GND. When connecting using an RS232C interface to a computer (DB9 socket (male type)) or USB-to-serial cable, use a male-to-female cross-serial cable. The wiring is as follows:

2-----2
3-----2
5-----5

Method for Setting and Options of the Electronic Scale:

Press the calibration button "CAL" and turn on the scale. After powering on, the scale will enter the setting mode. Press the "CAL" button to change the parameter that needs to be set (C1\C2\C3\C4\C5), press the "TARE" button to change the parameter value (0\1\2\3\4\5\6).

1. C1: Sensitivity Setting Options: 0, 1, 2, 3, and 4. The higher the value, the lower the sensitivity. Lower sensitivity provides better stability. The factory setting is 1 or 2.
2. C2: Filter Strength Setting Options: 0, 1, 2, and 3. The higher the value, the slower the response time. Slower response time improves stability. The factory setting is 1 or 2.
3. C3: Transmission Speed Setting (Factory setting is 2):
 - » 0 (9600 automatic continuous output)
 - » 1 (9600 automatic output after stabilization)
 - » 2 (600)
 - » 3 (1200)
 - » 4 (2400)
 - » 5 (4800)
 - » 6 (9600)
4. C4: Signal Transmission Setting
Press the "TARE" button to decrease the value, and press the "COUNT" button to increase the value. These data are the first data received during communication. The factory setting is 27.
5. C5: Function of the Unit Conversion Button (Factory setting is 0):
 - » 0 (unit conversion)
 - » 1 (print)

6. After setting the above parameters, press the "CAL" button again, and the scale will recheck the zero point.

7. Press the "TARE" button and turn the scale on again. After turning on the scale, all parameters and calibration data will be restored to the factory settings. The scale needs to be recalibrated before use.

5. PRECAUTIONS

1. The scale should not be used outdoors; it is designed for use at altitudes between 0 and 2000 meters, with an operating temperature range of 0 to 40°C, a relative humidity of $\leq 80\%$, and voltage fluctuations less than $\pm 10\%$.
2. The electronic scale is a precision instrument. Handle the object carefully when weighing, and avoid exceeding the maximum weighing range of the scale. Any form of overloading or impact can cause permanent damage to the scale, even if it is turned off. This also applies during use.
3. The working environment should be free of vibrations, interference, and corrosive substances.
4. Allow sufficient time for the scale to warm up before use (at least 15 minutes).
5. Use a dry cotton cloth for cleaning.

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INSIZE is a global manufacturer of measuring technology bringing technological innovations with representation in 75 countries around the world. Measuring instruments of the INSIZE brand represent the optimal solution without compromises to meet even the most demanding measuring needs.

INSIZE measuring instruments will pleasantly surprise you with:
| high-quality craftsmanship | reliable performance | affordable prices

INSIZE provides a complete portfolio of measuring instruments » calipers, height gauges, dial indicators, micrometers, roughness testers, hardness testers, measuring microscopes, optical measuring devices, profile projectors, testing machines, contour gauges, roundness measuring machines, thickness gauges, gages, protractors, force gauges, meters, scales, videoscopes, torque wrenches and accessories including stands for measuring instruments.

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